

Ventura 3-Wheel Scooter

HCPCS:



FIGURE A




FIGURE B




FIGURE C




FIGURE D



FIGURE E



FIGURE F

The Ventura 3 Wheel Scooter by Drive Medical comes with a Captain's Seat for superior comfort.

- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue (Figure A)
- 33 AH batteries for extended range
- Floorboard with 23" of foot space
- Ergonomic throttle control
- Height-adjustable swivel seat features fold-down backrest and adjustable padded armrests
- Seat position can be adjusted forward or backward
- Rear suspension and Anti-tip wheels
- Programmable controller
- Seat back reclines and features adjustable headrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust Delta tiller
- Headlight, Turn signals/hazard lights and Rearview mirrors
- Easy free wheel operation

PRODUCTS	Description
VENTURA320CS	VENTURA 3-WHL DELUXE 20"CAPT.SEAT,RED/BL
SPECIFICATIONS	Description
Warranty	24 Month Limited Warranty on Electronics
Product Weight	187 lb
Product Dimensions	48" (L) x 25" (W) x 45" (H)
Seat Dimensions	18" (D) x 20" (W)
Weight Capacity	400 lb
Brakes	Electromagnetic
Top Speed	6 mph
Batteries (pair)	12V x 33AH
Battery Charger	4A Offboard
Battery Weight (pair)	46 lbs.
Charger	4A Offboard
Climbing Angle	7.5°
Color	Red & Blue
Controller	PG 70A

SPECIFICATIONS	Description
Dimensions	48" (L) x 25" (W)
Ground Clearance	2.5"
HCPCS	N/A
Heaviest Piece Weight	60 lbs.
Maximum Range	16 miles
Motor	24V x 350W x 5100 RPM
Seat Type	Captain - Reclining High
Turning Radius	54"
Warranty on Batteries	12 Month Warranty on Batteries
Warranty on Frame	Limited Lifetime Warranty on Frame
Wheel Size	Front: 10" x 3" Rear: 10" x 3"

ACCESSORIES	Description
SF8030	Crutch Holder Set, SType, 1/ea, EA

PARTS	Description
VENT-74	Seat Plate, Capt. Seat Ventura 3&4, 1/ea, EA
VENT-73	Seat Plate, Folding Seat Ventura 3&4, EA
SEC2010N	HeadRestNew18"SunfireEC,1 ea, EA
VENT-65	Seat Slide Set, Ventura 3 & 4, 2/pkg, PKG
GT8072W32	Hazard Button,GladGT,Ventura,1/ea, EA
VENT-53	Mirror,Right, Ventura,Panther,1/ea, EA
VENT-54	Mirror,Left,Ventura,Panther,1/ea, EA
GT8072W35	FrontSignalLeftGladGT,Ventura,1/ea, EA
GT8072W36	FrontSignalRightGladGT,Ventura,1/ea, EA
GT8072W53	DirectionBulbGladGT,Ventura,1/ea, EA
VENT-55	Basket Holder, Ventura, 1/ea, EA
VENT-61	Front Decoration, Red, Ventura3&4, 1/ea, EA
VENT-62	Front Decoration, Blue, Ventura3&4, 1/ea, EA
VENT-46	Coupler, Ventura 3, 1/ea, EA
VENT-49	FrontWheelAssemblyVentura3,1/ea, EA
VENT-45	FrontShroudVentura3,1/ea, EA
VENT-50	FrontWheelAxle,Ventura 3,1/ea, EA
VENT-48	FrontForkAssemblyVentura3,1/ea, EA
VENT-58	FrontWheelW/PneumaticTire,Ventura3,1/ea, EA
VENT-59	OuterTireFront&RearWheelPneumTireVent3&4, EA