

Model	Maxima <sup>1</sup>
Part Number	3W: SC901 4W: SC941
Weight Capacity	500 lbs.
Maximum Speed <sup>2</sup>	Up to 5.8 mph
Ground Clearance <sup>3</sup>	3W: 2.5" at motor; 4W: 2.25" at motor
Turning Radius <sup>3</sup>	3W: 49"; 4W: 53"
Overall Length <sup>3</sup>	47.5"
Overall Width <sup>3</sup>	27" (not including armrest)
Seat-to-Ground Height <sup>3</sup>	(3W/4W): 22.1" – 24.1" 22.32" – 28.57" with Power Elevating Seat
Seat-to-Deck Height <sup>3</sup>	(3W/4W): 17.13" – 19.13" 17.36" – 23.61" with Power Elevating Seat
Front Tires <sup>3</sup>	3W: 10.5" x 3.6" solid 4W: 10" x 3.6" solid
Rear Tires <sup>3</sup>	10.75" x 3.6" solid
Suspension	N/A
Range Per Charge <sup>2,5</sup>	3W: Up to 19.8 miles 4W: Up to 19.5 miles
Total Weight With and Without Batteries <sup>7</sup>	3W (w/out batteries): 183.5 lbs. 3W (w/batteries): 256.5 lbs. 4W (w/out batteries): 196.8 lbs. 4W (w/batteries): 269.8 lbs.
Heaviest Piece When Disassembled <sup>3</sup>	N/A (unit does not disassemble)
Standard Seating <sup>3</sup>	Type: Depth-adjustable, high-back Weight: 53.2 lbs. Material: Black vinyl Dimensions: 22"x 18"-20" (depth adjustable) Back height: 19.75"
Drive System	Sealed transaxle, 24-volt DC motor
Dual Braking System	Regenerative and electromechanical
Battery Requirements <sup>4,6</sup>	Size: (2) Gel, 12V, 58Ah, NF-22 Weight: 36.5 lbs. each
Battery Charger	Off-board, 5-amp
Warranty	Frame: lifetime limited Drivetrain: 2-year limited Electronics: 2-year limited Batteries: 6-month limited
Color Option	Garnet Red (matte) Ocean Blue (matte)

3W: HCPCS Code: K0802 (Group 1, Very Heavy Duty)  
4W: HCPCS Code: K0808 (Group 2, Very Heavy Duty)

- 1 Pride® FDA Class II Medical Devices are designed to aid individuals with mobility impairments.
- 2 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition and tire condition. These specifications can be subject to a variance of (+/- 10%).
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or - 3%).
- 4 AGM or gel cell type required.
- 5 Tested in accordance with ANSI/RESNA, WC Vol2, section 4 & ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- 6 Battery weight may vary based on manufacturer.
- 7 Includes standard seat weight.



(US) 800-800-8586 • (Canada) 888-570-1113

pridemobility.com |    

HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. DME providers are responsible for determining appropriate billing codes when submitting for insurance reimbursement.

The information contained herein is correct at time of publication; we reserve the right to alter specifications without prior notice.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Special order seats or non-standard dimensions will not be fold-flat and depth-adjustable.

MAXIMA\_071421

# MAXIMA<sup>1</sup>

## ULTRA HEAVY-DUTY FOR MAXIMUM PERFORMANCE



FDA Class II  
Medical Device<sup>1</sup>

LIVE YOUR BEST®

pridemobility.com |    

# HEAVY-DUTY PERFORMANCE

22" x 18"-20"  
High back black vinyl, depth adjustable, limited recline and swivel seat

Rear taillights

Black, flat-free, non-scuffing tires

Bumpers

Delta tiller

Tiller lever

USB charger

Lower deck-to-ground

Rearview mirror

Standard front basket

Full LED light package

4W — FDA Class II Medical Device<sup>1</sup>



Full LED light package

Infinite position tiller lever (mech lock)



3W — FDA Class II Medical Device<sup>1</sup>

## ALL NEW REDESIGN!

### NEW! Tiller adjustment

- Infinite position tiller lever (mech lock)

### NEW! Seat

- 22"x 18"-20" High back depth adjustable
- 106° limited recline
- 4-point swivel (no sliders)
- 7" headrest with 3" of adjustment

**NEW! USB charger** built into the tiller for convenient charging of smart phones and portable electronic devices

**NEW! Full LED light package** includes front headlights, turn signals and rear taillights

**NEW! Ambient voltmeter** screen brightness automatically adjusts based on light intensity


**NEW! Delta tiller** with ergonomic wraparound handles lets you operate the scooter with one hand

- Longer per-charge range
- Standard with front and side bumper (3W) and side bumper (4W)
- XLR battery charger on tiller

### Console

- Ambient voltmeter
- On/off key switch
- Horn
- Speed control
- Light switch
- Turn signal button

### Colors



Garnet Red (matte)

Ocean Blue (matte)

### Option

- Power elevating seat available on 3 & 4 wheel
- 6.25" of additional height to access cabinets, countertops and other out-of-reach items

### Accessories

- Additional rearview mirror
- Cup holder (B)
- Forearm crutch holder
- Front bumper (4w)
- Oxygen holder
- Quad cane holder
- Rear basket
- Rear bumper (3w, 4w)
- Saddle bag
- Safety flag
- Single/double cane/crutch holder
- Weather covers
- Wishbone crutch holder

### Military Patches

