

# Enduro XL3

Model Number: 778XLSBN

Medicare Billing Code: K0899




The Enduro XL3 is a heavy duty three-wheel scooter with a 500 lb. weight capacity. It has a delta style tiller which houses the controls for the lighting system, turn signals, wig wag and horn. It has a dual purpose wig wag which allows the user to push forward or pull back on the drive lever depending on dexterity for ease of driving. The unit comes standard with a 20in. seat that can accommodate most riders, and has swivel and recline features for added comfort and ease of use.

- + 500lb Weight Capacity
- + Slide & Swivel Seating
- + Full Suspension System
- + Full Lighting Package
- + Adjustable Tiller



Available Color:



## SPECIFICATIONS

WEIGHT CAPACITY	500 lbs
WIDTH	27.5"
LENGTH	55"
UNIT WEIGHT <sup>1</sup>	255 lbs
MAX SPEED <sup>2</sup>	9 mph

GROUND CLEARANCE	4"
SAFE CLIMBING ANGLE	8°   14%
FRONT WHEELS	12.5" Pneumatic
REAR WHEELS	14" Pneumatic
SUSPENSION TYPE	Full

BATTERY OPTIONS	50Ah   60Ah   75Ah
BATTERY WEIGHTS	50 lbs ea
RANGE <sup>3</sup> UP TO	24   29   37 miles
CHARGER	8A Off-board
TILLER TYPE	Delta

Note: <sup>1</sup> Maximum unit weight dependent on seating configuration and battery choice, standard equipment weight shown. <sup>2</sup> Maximum speed may vary, user weight and road grade may affect top speed. <sup>3</sup> Range can vary depending on level of charge, user weight and road conditions. HCPCS codes provided should not be considered as legal advice, and does not guarantee reimbursement.